



# MYB75 Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA863725LD01DOA
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9FE25
<b>Immunogen</b>	Recombinant Arabidopsis thaliana Transcription factor MYB75 protein (91-167AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Arabidopsis thaliana
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Transcription factor MYB75 (Myb-related protein 75) (AtMYB75) (Production of anthocyanin pigment 1 protein) (Suc-induced anthocyanin accumulation locus 1), MYB75, PAP1 SIAA1
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Research Area</b>	Others
<b>Target Names</b>	MYB75
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.