



TBC1D23 Antibody, FITC conjugated

Product Code	CSB-PA865143LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NUY8
Immunogen	Recombinant Human TBC1 domain family member 23 protein (510-605AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	TBC1 domain family member 23 (HCV non-structural protein 4A-transactivated protein 1), TBC1D23, NS4ATP1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TBC1D23
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.