



CLN8 Antibody, HRP conjugated

Product Code	CSB-PA866210LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBY8
Immunogen	Recombinant Human Protein CLN8 protein (247-286AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein CLN8, CLN8, C8orf61
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	CLN8