



RAI14 Antibody, HRP conjugated

Product Code	CSB-PA868376LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P0K7
Immunogen	Recombinant Human Ankyrin protein (242-346AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ankyrin (Ankyrin repeat and coiled-coil structure-containing protein) (Novel retinal pigment epithelial cell protein) (Retinoic acid-induced protein 14), RAI14, KIAA1334 NORPEG
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RAI14
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.