



RAI14 Antibody, HRP conjugated

Product Code	CSB-PA868376LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P0K7
Immunogen	Recombinant Human Ankycorbin protein (242-346AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ankycorbin (Ankyrin repeat and coiled-coil structure-containing protein) (Novel retinal pigment epithelial cell protein) (Retinoic acid-induced protein 14), RAI14, KIAA1334 NORPEG
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RAI14