



MIEN1 Antibody, Biotin conjugated

Product Code	CSB-PA874795LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BRT3
Immunogen	Recombinant Human Migration and invasion enhancer 1 protein (2-112AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Migration and invasion enhancer 1 (HBV X-transactivated gene 4 protein) (HBV XAg-transactivated protein 4) (Protein C35), MIEN1, C17orf37 RDX12 XTP4
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	MIEN1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.