



MSL2 Antibody, Biotin conjugated

Product Code	CSB-PA875722LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HCI7
Immunogen	Recombinant Human E3 ubiquitin-protein ligase MSL2 protein (341-473AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase MSL2 (EC 2.3.2.-) (E3 ubiquitin-protein transferase MSL2) (Male-specific lethal 2-like 1) (MSL2-like 1) (Male-specific lethal-2 homolog) (MSL-2) (Male-specific lethal-2 homolog 1) (RING finger protein 184), MSL2, KIAA1585 MSL2L1 RNF184
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	MSL2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.