



KIF13A Antibody, HRP conjugated

Product Code	CSB-PA880944LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H1H9
Immunogen	Recombinant Human Kinesin-like protein KIF13A protein (1383-1537AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kinesin-like protein KIF13A (Kinesin-like protein RBKIN), KIF13A, RBKIN
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KIF13A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.