



GPR27 Antibody, HRP conjugated

Product Code	CSB-PA885717LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NS67
Immunogen	Recombinant Human Probable G-protein coupled receptor 27 protein (203-285AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable G-protein coupled receptor 27 (Super conserved receptor expressed in brain 1), GPR27, SREB1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GPR27