



GPR27 Antibody, HRP conjugated

Product Code	CSB-PA885717LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NS67
Immunogen	Recombinant Human Probable G-protein coupled receptor 27 protein (203-285AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable G-protein coupled receptor 27 (Super conserved receptor expressed in brain 1), GPR27, SREB1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GPR27
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.