



HHLA2 Antibody, HRP conjugated

Product Code	CSB-PA891555LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UM44
Immunogen	Recombinant Human HERV-H LTR-associating protein 2 protein (172-258AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	HERV-H LTR-associating protein 2 (Human endogenous retrovirus-H long terminal repeat-associating protein 2), HHLA2
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	HHLA2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.