



HLA2 Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA891555LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9UM44 |
| Immunogen | Recombinant Human HERV-H LTR-associating protein 2 protein (172-258AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | HERV-H LTR-associating protein 2 (Human endogenous retrovirus-H long terminal repeat-associating protein 2), HHLA2 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Immunology |
| Target Names | HLA2 |