



FHOD1 Antibody, HRP conjugated

Product Code	CSB-PA896540LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y613
Immunogen	Recombinant Human FH1/FH2 domain-containing protein 1 protein (338-437AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	FH1/FH2 domain-containing protein 1 (Formin homolog overexpressed in spleen 1) (FHOS) (Formin homology 2 domain-containing protein 1), FHOD1, FHOS FHOS1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FHOD1