




IFT52 Antibody

Product Code	CSB-PA896883LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y366
Immunogen	Recombinant Human Intraflagellar transport protein 52 homolog protein (325-434AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Intraflagellar transport protein 52 homolog (Protein NGD5 homolog), IFT52, C20orf9 NGD5
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	IFT52

Image



Western Blot

Positive WB detected in: A549 whole cell lysate, HepG2 whole cell lysate, PC-3 whole cell lysate
All lanes: IFT52 antibody at 5.1µg/ml
Secondary
Goat polyclonal to rabbit IgG at 1/50000 dilution
Predicted band size: 50 kDa
Observed band size: 50 kDa