




IFT52 Antibody

| | |
|----------------------------|-----------------------------------------------------------------------------------------|
| Product Code | CSB-PA896883LA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9Y366 |
| Immunogen | Recombinant Human Intraflagellar transport protein 52 homolog protein (325-434AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA, WB; Recommended dilution: WB:1:500-1:5000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Intraflagellar transport protein 52 homolog (Protein NGD5 homolog), IFT52, C20orf9 NGD5 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Others |
| Target Names | IFT52 |

Image



Western Blot

Positive WB detected in: A549 whole cell lysate, HepG2 whole cell lysate, PC-3 whole cell lysate

All lanes: IFT52 antibody at 5.1µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 50 kDa

Observed band size: 50 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.