



IFT52 Antibody, FITC conjugated

Product Code	CSB-PA896883LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y366
Immunogen	Recombinant Human Intraflagellar transport protein 52 homolog protein (325-434AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Intraflagellar transport protein 52 homolog (Protein NGD5 homolog), IFT52, C20orf9 NGD5
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	IFT52
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.