



# IFT52 Antibody, Biotin conjugated

Product Code	CSB-PA896883LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y366
Immunogen	Recombinant Human Intraflagellar transport protein 52 homolog protein (325-434AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Intraflagellar transport protein 52 homolog (Protein NGD5 homolog), IFT52, C20orf9 NGD5
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	IFT52