



# Goat Anti-Rabbit IgG(H+L) Antibody; HRP conjugated

<b>Product Code</b>	CSB-PA564648
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Raised In</b>	Goat
<b>Species Reactivity</b>	Rabbit
<b>Tested Applications</b>	ELISA, WB, IHC, ICC; Recommended dilution: ELISA:1:500-1:10000, WB:1:3000-1:5000, IHC:1:50-1:500, ICC:1:50-1:500
<b>Form</b>	Liquid
<b>Conjugate</b>	HRP
<b>Storage Buffer</b>	Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
<b>Isotype</b>	IgG
<b>Product Type</b>	Secondary Antibody
<b>Immunogen Species</b>	Oryctolagus cuniculus (Rabbit)