



Goat Anti-Rabbit IgG(H+L) Antibody; HRP conjugated

Product Code	CSB-PA564648
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Raised In	Goat
Species Reactivity	Rabbit
Tested Applications	ELISA, WB, IHC, ICC; Recommended dilution: ELISA:1:500-1:10000, WB:1:3000-1:5000, IHC:1:50-1:500, ICC:1:50-1:500
Form	Liquid
Conjugate	HRP
Storage Buffer	Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Isotype	IgG
Product Type	Secondary Antibody
Immunogen Species	Oryctolagus cuniculus (Rabbit)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.