



# Recombinant Human Neurturin protein (NRTN)

<b>Product Code</b>	CSB-AP002961HU
<b>Uniprot No.</b>	Q99748
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>96% as determined by SDS-PAGE and HPLC.
<b>Sequence</b>	ARLGARPCGL RELEVRVSEL GLGYASDETV LFRYCAGACE AAARVYDLGL RRLRQRRRLR RERVRAQPCC RPTAYEDEVS FLDAHSRYHT VHEL SARECA CV
<b>Source</b>	E.coli
<b>Target Names</b>	NRTN
<b>Tag Info</b>	N-terminal His-tagged Tag-Free
<b>Mol. Weight</b>	23.4 kDa
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.