



Recombinant Human Metalloproteinase inhibitor 2 protein (TIMP2)

Product Code	CSB-AP000161HU
Uniprot No.	P16035
Alias	CSC-21K, TIMP-2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% as determined by SDS-PAGE and HPLC.
Sequence	CSCSPVHPQQ AFCNADVIR AKAVSEKEVD SGNDIYGNPI KRIQYEIKQI KMFKGPEKDI EFIYTAPSSA VCGVSLDVGG KKEYLIAGKA EGDGKMHITL CDFIVPWDTL STTQKKSLNH RYQMGCECKI TRCPMIPCYI SSPDECLWMD WVTEKNINGH QAKFFACIKR SDGSCAWYRG AAPPKQEFLD IEDP
Source	E.coli
Target Names	TIMP2
Tag Info	N-terminal His-tagged Tag-Free
Mol. Weight	21.8 kDa
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.