



# Recombinant Human Metalloproteinase inhibitor 2 protein (TIMP2)

<b>Product Code</b>	CSB-AP000161HU
<b>Uniprot No.</b>	P16035
<b>Alias</b>	CSC-21K, TIMP-2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% as determined by SDS-PAGE and HPLC.
<b>Sequence</b>	CSCSPVHPQQ AFCNADVIR AKAVSEKEVD SGNDIYGNPI KRIQYEIKQI KMFKGPEKDI EFIYTAPSSA VCGVSLDVGG KKEYLIAGKA EGDGKMHTL CDFIVPWDTL STTQKKSLSNH RYQMGCECKI TRCPMIPCYI SSPDECLWMD WVTEKNINGH QAKFFACIKR SDGSCAWYRG AAPPKQEFLD IEDP
<b>Source</b>	E.coli
<b>Target Names</b>	TIMP2
<b>Tag Info</b>	N-terminal His-tagged Tag-Free
<b>Mol. Weight</b>	21.8 kDa
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.