



Recombinant *Vibrio proteolyticus* protein (AMPX)

Product Code	CSB-AP000191VCN
Uniprot No.	Q01693
Product Type	Recombinant Protein
Immunogen Species	<i>Vibrio proteolyticus</i>
Purity	>98% as determined by SDS-PAGE and HPLC.
Sequence	MPPITQQATV TAWLPQVDAS QITGTISSLE SFTNRFYTTT SGAQASDWIA SEWQALSASL PNASVKQVSH SGYNQKSVVM TITGSEAPDE WIVIGGHLDS TIGSHTNEQS VAPGADDAS GIAAVTEVIR VLSENNFQPK RSIAFMAYAA EEVGLRGSQD LANQYKSEGK NVVSALQLDM TNYKGS AQDV VFITDYTDSN FTQYLTQLMD EYLPSLTYGF DTCGYACSDH ASWHNAGYPA AMPFESKFND YNPRIHTTQD TLANS DPTGS HAKKFTQLGL AYA IEMGSAT G
Source	E.coli
Tag Info	N-terminal His-tagged Tag-Free
Mol. Weight	31.4 kDa
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.