



Recombinant Heloderma suspectum Exendin-4

Product Code	CSB-AP000231HVJ
Uniprot No.	P26349
Product Type	Recombinant Protein
Immunogen Species	Heloderma suspectum?
Purity	≥96% as determined by SDS-PAGE and HPLC.
Sequence	HGEGTFTSDL SKQMEEEEAVR LFIEWLKNNGG PSSGAPPPS
Source	E.coli
Tag Info	N-terminal His-tagged Tag-Free
Mol. Weight	4.2 kDa
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.