



Recombinant Human MHC class I polypeptide-related sequence A protein (MICA)

Product Code	CSB-AP000261HU
Uniprot No.	Q29983
Alias	MIC-A,
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% as determined by SDS-PAGE and HPLC.
Sequence	EPHSLRYNLT VLSWDGSVQS GFLAEVHLDG QPFLRYDRQK CRAKPQGQWA EDVLGNKTWD RETRDLTGNG KDLRMTLAHI KDQKEGLHSL QEIRVCEIHE DNSTRSSQHF YYDGELFLSQ NLETEEWTV P QSSRAQTLAM NVRNFLKEDA MKTKTHYHAM HADCLQELRR YLESGVVLRR TVPPMVNVTR SEASEGNITV TCRASSFYPR NIILTW RQDG VLSHDTQQW GDVLPDGN GT YQTWVATRIC RGEEQRFTCY MEHSGNHSTH PVPSGKVLVL QSHW
Source	E.coli
Target Names	MICA
Tag Info	N-terminal His-tagged Tag-Free
Mol. Weight	32.8 kDa
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.