



Recombinant Human MHC class I polypeptide-related sequence B protein (MICB)

Product Code	CSB-AP000271HU
Uniprot No.	Q29980
Alias	MIC-B,
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% as determined by SDS-PAGE and HPLC.
Sequence	AEPHSLRYNL MVLSQDESVQ SGFLAEGHLD GQPFLRYDRQ KRRAKPQGQW AEDVLGAKTW DTETEDLTEN GQDLRRTLTH IKDQKGGGLHS LQEIRVCEIH EDSSTRGSRH FYYDGELFLS QNLETQESTV PQSSRAQTLA MNVTNFWKED AMKTKTHYRA MQADCLQKLQ RYLKSGVAIR RTVPPMVNVT CSEVSEGNIT VTCRASSFYP RNITLTWRQD GVSLSHNTQQ WGDVLPDGNG TYQTWVATRI RQGEEQRFTC YMEHSGNHGT HPVPSGKVLV LQSQRTD
Source	E.coli
Target Names	MICB
Tag Info	N-terminal His-tagged Tag-Free
Mol. Weight	32.8 kDa
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.