



Recombinant Human Nucleobindin-2 protein (NUCB2)

Product Code	CSB-AP000391HU
Uniprot No.	P80303
Alias	DNA-binding protein NEFA, Prepronesfatin,
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% as determined by SDS-PAGE and HPLC.
Sequence	VPIDIDKTKV QNIHPVESAK IEPDPTGLYY DEYLKQVIDV LETDKHFREK LQKADIEEIK SGRLSKELDL VSHHVRTKLD EL
Source	E.coli
Target Names	NUCB2
Tag Info	N-terminal His-tagged Tag-Free
Mol. Weight	9.6 kDa
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.