




Recombinant Human Receptor-type tyrosine-protein phosphatase zeta (PTPRZ1), partial

Product Code	CSB-EP019068HUa0
Relevance	Protein tyrosine phosphatase that negatively regulates oligodendrocyte precursor proliferation in the embryonic spinal cord. Required for normal differentiation of the precursor cells into mature, fully myelinating oligodendrocytes. May play a role in protecting oligodendrocytes against apoptosis. May play a role in the establishment of contextual memory, probably via the dephosphorylation of proteins that are part of important signaling cascades
Abbreviation	Recombinant Human PTPRZ1 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P23471
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	IGWSYTGALNQKNWGGKYPCTCNPKQSPINIDEDLTQVNVNLKKLKFQGWDK TSLENTFIHNTGKTVEINLTNDYRVSGGVSEMVFKASKITFHWGKCNMSSDGS EHSLEGQKFPLEMQIYCFDADRFSSFEEAVKGKGKLRALSILFEVGTEENLDFK AIIDGVESVSRFGKQAALDPFILLNLLPNSTDKYYIYNGSLTSPPCDTVDWIVF KDTVSISESQ LAVFCEVLTMQQSGYVMLMDY LQNNFREQQYKFSRQVFSSY
Research Area	Neuroscience
Source	E.coli
Target Names	PTPRZ1
Protein Names	Protein-tyrosine phosphatase receptor type Z polypeptide 1 (Protein-tyrosine phosphatase receptor type Z polypeptide 2) (R-PTP-zeta-2) (HTPZP2) (PTPRZ) (PTPRZ2) (PTPZ)
Expression Region	36-300aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	34.1 kDa
Protein Length	partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.