



rbtD Antibody, Biotin conjugated

Product Code	CSB-PA360400HD01EKN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00335
Immunogen	Recombinant Enterobacter aerogenes Ribitol 2-dehydrogenase protein (1-249AA)
Raised In	Rabbit
Species Reactivity	Enterobacter aerogenes
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ribitol 2-dehydrogenase (RDH) (EC 1.1.1.56), rbtD
Immunogen Species	Enterobacter aerogenes
Research Area	Others
Target Names	rbtD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.