




Recombinant Human Heterogeneous nuclear ribonucleoproteins A2/B1 (HNRNPA2B1)

Product Code	CSB-EP010602HUb0
Abbreviation	Recombinant Human HNRNPA2B1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	I6L957
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MEREKEQFRKLFIGGLSFETTEESLRNYYEQWGKLTDCVVMRDPASKRSRGF GFVTFSSMAEVDAAAMAARPHSIDGRVVEPKRAVAREESGKPGAHVTVKKLFFV GGIKEDTEEHHLRDYFEEYGKIDTIEIITDRQSGKKRGFGFVTFDDHDPVDKIVL QKYHTINGHNAEVRKALSRQEMQEDLEVAILEVAPVMEEEEEEDMVVEDLDMA TRVGATEVVMTTMEEEIMEVEITMILEIITSNLLTTVQ
Research Area	Epigenetics and Nuclear Signaling
Source	E.coli
Target Names	HNRNPA2B1
Expression Region	1-249aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Mol. Weight	31.9 kDa
Protein Length	Full Length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.