





Recombinant Human Desmoglein-3 (DSG3), partial


Product Code	CSB-EP007205HUa0
Abbreviation	Recombinant Human DSG3 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P32926
Storage Buffer	Tris-based buffer, 50% glycerol
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	EWVKFAKPCREGEDNSKRNP IAKITSYQATQKITYRISGV GIDQPPFGIFVVDKNTG DINITAIVDREETPSFLIT CRALNAQGLDVEKPLIL TVKILDINDNPPVFSQQ IFMGEIEENSASNSLVM ILNATDADEPNHLNSKI AFKIVSQEPAGTPMFLL SRNTGEVRTLTNSLDRE QASSYRLVSGADKDGEG LSTQCECNKVKDVNDNF PMFRDSQYSARIEENIL SSELLRFQVTDLDEEYTD NWLAVYFFTSGNEGNWF EIQTDPRTNEGILKVVK ALDYEQLQSVKLSIAVK NKAEFHQSVISRYRVQS TPVTIQVINVREGIAFR PASKTFTVQKGISSKKLV DYILGTYQAIDEDTNKA ASNKYYVMGRNDGGYLM IDSKTAEIKFVKNMNRD STFIVNKTITAEVLAIDE YTGKTSTGTVYVRVPDF NDNCPTAVLEKDAVCSS SPSVVVSARTLNNRYTG PYTFALEDQPVKLPAV WSITTLNATSALLRAQE QIPPGVYHISLVLTDSQ NNRCEMPRSLTLEVCQ CDNRGICGTSYPTTSPG TRYGRPHSGR
Research Area	Others
Source	E.coli
Target Names	DSG3
Expression Region	50-615aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	67.0 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Based on the SEQUEST from database of E.coli host and target protein, the LC-MS/MS Analysis result of CSB-EP007205HUa0 could indicate that this peptide derived from E.coli-expressed Homo sapiens (Human) DSG3.



Based on the SEQUEST from database of E.coli host and target protein, the LC-MS/MS Analysis result of CSB-EP007205HUa0 could indicate that this peptide derived from E.coli-expressed Homo sapiens (Human) DSG3.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.