




Recombinant *Saccharomyces cerevisiae* Acetyl-CoA carboxylase (ACC1), partial

Product Code	CSB-EP2725SVG1
Abbreviation	Recombinant <i>Saccharomyces cerevisiae</i> ACC1 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q00955
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	TGLVSVDDDIYQKGCCTSPEDGLQKAKRIGFPVMIKASEGGGGKGIRQVEREE DFIALYHQAANEIPGSPIFIMKLAGRARHLEVQLLADQYGTNISLFGRDCSVQR RHQKIIIEAPVTIAKAETFHEMEKAAVRLGKLVGYVSAGTVEYLYSHDDGKFYF LELNPRQLQVEHPTTEMVSGVNLPAALQIAMGIPMHRISDIRTLYGMNPHSASE IDFEFKTQDATKKQRRPIPKGHCTACRITSEDPNDGFKPSGGTLHELNFRSSS NVWGYFSVGNNNGNIHSFSDSQFGHIFAFGENRQASRKHMVVALKELSIRGDF RTTVEYLIKLLLETEDFEDNTITTGWLDDLITHKMTAEKPDPTLAVICGAATKAFLA SEEARHKYIESLQKGQVLSKDLLQTMFPVDFIHEGKRYKFTVAKSGNDRYTLFI NGSKCDIILRQLSDGGLLIAIGGKSHTIYWKEEVAATRLSVDSMTTLLEVENDPT QLRTPSPGKLVKFLVENGEHIIKGQPYAEIEVMKMQMPLVVSQENGIVQLLKQP GSTIVAGDIMAIMTLDDPSKV
Research Area	Cancer
Source	E.coli
Target Names	ACC1
Expression Region	217-775aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Mol. Weight	67.6 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.