



Recombinant Human Peroxisome proliferator-activated receptor alpha (PPARA)

Product Code	CSB-BP018421HU
Abbreviation	Recombinant Human PPARA protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q07869
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	MVDTESPLCPLSPLEAGDLESPLSEEFLLQEMGNIQEISQSIGEDSSGSFGFTEY QYLGSCPGSDGSVITDTLSPASSPSSVTYPVVPGSVDESPSGALNIECRICGD KASGYHYGVHACEGCKGFFRRTIRLKLVDKCDRSCIKKKNRNCQYCRFH KCLSVGMASHNAIRFGRMPRSEKAKLKAIEILTCEHDIEDSETADLKSLAKRIYEA YLKNFNMNKVKARVILSGKASNNPPFVIHDMETLCMAEKTLVAKLVANGIQNKE AEVRIFHCCQCTSVETVTELTEFAKAIPGFANLDLNDQVTLTKYGVYEAIFAMLS SVMNKDGMVLVAYGNFITREFLKSRLKPFCDIMEPKFDFAMKFNALDDSDIS LFVAAIICCGDRPGLLNVDGHIEMQEGIVHVLRLHLQSNHPDDIFLFPKLLQKMA DLRQLVTEHAQLVQIIKKTESDAALHPLLQEYRDMY
Research Area	Cancer
Source	Baculovirus
Target Names	PPARA
Expression Region	1-468aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	56.2 kDa
Protein Length	Full Length
Image	

