



Recombinant Human PiggyBac transposable element-derived protein 5 (PGBD5)

Product Code	CSB-BP854034HU
Abbreviation	PGBD5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N414
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAEGGGGARRRAPALLEAARARYESLHISDDVFGESGPDSSGGNPFYSTSAAS RSSSAASSDDEREPPGPPGAAPPPRAPDAQEPEEDEAGAGWSAALRDRPP PRFEDTGGPTRKMPPSASAVDFQFLVDPNVLKNMVVQTNMYAKKFQERFG SDGAWVEVTLTEMKAFLGYMISTSISHCESVLSIWGGGFYSNRSLALVMSQAR FEKILKYFHVVAFRSSQTTHGLYKVPFLDSLQNSFDSAFRPSQTQVLHEPLID EDPVFIATCTERELRKRKRKFSWLWRQCSSTGFIIQIYVHLKEGGGPDGLDAL KNKPQLHSMVARSLCRNAAGKNYIIFTGPSITSLTLFEEFEKQGIYCCGLLRARK SDCTGLPLSMLTNPATPPARGQYQIKMKGNMSLICWYNKGHFRFLTNAYSVP QQGVIIKRKSGEIPCPLAVEAFAAHLASYICRYDDKYSKYFISHKPNKTWQQVFW FAISAIINNAYILYKMSDAYHVKRYRAQFGERLVRELLGLEDasPTH
Source	Baculovirus
Target Names	PGBD5
Protein Names	Recommended name: PiggyBac transposable element-derived protein 5
Expression Region	1-524
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.