



BPIFB3 Antibody, HRP conjugated

Product Code	CSB-PA003693LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P59826
Immunogen	Recombinant Human BPI fold-containing family B member 3 protein (248-352AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	BPI fold-containing family B member 3, Ligand-binding protein RYA3, Long palate, lung and nasal epithelium carcinoma-associated protein 3, BPIFB3, C20orf185, LPLUNC3
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	BPIFB3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.