



CNGA1 Antibody, Biotin conjugated

Product Code	CSB-PA005642LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29973
Immunogen	Recombinant Human cGMP-gated cation channel alpha-1 protein (586-690AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	cGMP-gated cation channel alpha-1, Cyclic nucleotide-gated cation channel 1, Cyclic nucleotide-gated channel alpha-1, CNG channel alpha-1, CNG-1, CNG1, Cyclic nucleotide-gated channel, photoreceptor, Rod photoreceptor cGMP-gated channel subunit alpha, CNGA1, CNCG, CNCG1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal transduction
Target Names	CNGA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.