



LIPC Antibody, FITC conjugated

Product Code	CSB-PA012974LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11150
Immunogen	Recombinant Human Hepatic triacylglycerol lipase protein (23-217AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hepatic triacylglycerol lipase, HL, Hepatic lipase, Lipase member C, LIPC, HTGL
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	LIPC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.