



LIPC Antibody, Biotin conjugated

Product Code	CSB-PA012974LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11150
Immunogen	Recombinant Human Hepatic triacylglycerol lipase protein (23-217AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hepatic triacylglycerol lipase, HL, Hepatic lipase, Lipase member C, LIPC, HTGL
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	LIPC