



ntpcr Antibody, Biotin conjugated

Product Code	CSB-PA24179D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B8K0C4
Immunogen	Recombinant Zebrafish Nucleoside-triphosphatase, cancer-related protein (1-185AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nucleoside-triphosphatase, cancer-related, ntpcr
Immunogen Species	Zebrafish
Target Names	ntpcr
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.