



Ita Antibody, FITC conjugated

Product Code	CSB-PA24269C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q4W897
Immunogen	Recombinant Zebrafish Lymphotoxin alpha (TNF superfamily, member 1) protein (116-203AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lymphotoxin alpha (TNF superfamily, member 1), Ita
Immunogen Species	Zebrafish
Target Names	Ita
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.