



Ita Antibody, FITC conjugated

Product Code	CSB-PA24269C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q4W897
Immunogen	Recombinant Zebrafish Lymphotoxin alpha (TNF superfamily, member 1) protein (116-203AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lymphotoxin alpha (TNF superfamily, member 1), Ita
Immunogen Species	Zebrafish
Target Names	Ita