



lats2 Antibody, HRP conjugated

Product Code	CSB-PA24329B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B0S8I4
Immunogen	Recombinant Zebrafish Large tumor suppressor kinase 2 protein (1022-1140AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Large tumor suppressor kinase 2, lats2
Immunogen Species	Zebrafish
Target Names	lats2