



lats2 Antibody, Biotin conjugated

Product Code	CSB-PA24329D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B0S8I4
Immunogen	Recombinant Zebrafish Large tumor suppressor kinase 2 protein (1022-1140AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Large tumor suppressor kinase 2, lats2
Immunogen Species	Zebrafish
Target Names	lats2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.