



# lats2 Antibody, Biotin conjugated

|                     |  |
|---------------------|--|
| Product Code        | CSB-PA24329D0Rb  |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                    |
| Uniprot No.         | B0S8I4   |
| Immunogen           | Recombinant Zebrafish Large tumor suppressor kinase 2 protein (1022-1140AA)      |
| Raised In           | Rabbit   |
| Species Reactivity  | Zebrafish  |
| Tested Applications | ELISA  |
| Form                | Liquid   |
| Conjugate           | Biotin   |
| Storage Buffer      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified   |
| Isotype             | IgG  |
| Clonality           | Polyclonal   |
| Alias               | Large tumor suppressor kinase 2, lats2   |
| Immunogen Species   | Zebrafish  |
| Target Names        | lats2  |