



notchl Antibody, HRP conjugated

Product Code	CSB-PA24369B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	X1WEZ2
Immunogen	Recombinant Zebrafish Notch homolog,-like protein (683-784AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Notch homolog,-like, notchl
Immunogen Species	Zebrafish
Target Names	notchl
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.