



egf Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA24449C0Rb |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | A0A140LGI0 |
| Immunogen | Recombinant Zebrafish Pro-epidermal growth factor protein (1063-1176AA) |
| Raised In | Rabbit |
| Species Reactivity | Zebrafish |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Pro-epidermal growth factor, egf |
| Immunogen Species | Zebrafish |
| Target Names | egf |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |