



met Antibody, Biotin conjugated

Product Code	CSB-PA24469D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	F1QLR7
Immunogen	Recombinant Zebrafish MET proto-oncogene, receptor tyrosine kinase protein (76-179AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	MET proto-oncogene, receptor tyrosine kinase, met
Immunogen Species	Zebrafish
Target Names	met
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.