




ptprsa Antibody

Product Code	CSB-PA24489A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	F1QLA4
Immunogen	Recombinant Zebrafish Protein tyrosine phosphatase, receptor type, s, a protein (269-424AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein tyrosine phosphatase, receptor type, s, a, ptpsa
Immunogen Species	Zebrafish
Target Names	ptpsa

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: ptpsa antibody at 3.1µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 48 kDa

Observed band size: 60 kDa