



ptprsa Antibody, Biotin conjugated

Product Code	CSB-PA24489D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	F1QLA4
Immunogen	Recombinant Zebrafish Protein tyrosine phosphatase, receptor type, s, a protein (269-424AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein tyrosine phosphatase, receptor type, s, a, ptprsa
Immunogen Species	Zebrafish
Target Names	ptprsa
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.