



olig2 Antibody, Biotin conjugated

Product Code	CSB-PA24549D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A8W1Q8
Immunogen	Recombinant Zebrafish Oligodendrocyte lineage transcription factor 2 protein (150-273AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Oligodendrocyte lineage transcription factor 2, olig2
Immunogen Species	Zebrafish
Target Names	olig2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.